

# Axen™ Total DNA mini Kit

## Product description

Axen™ Total DNA mini Kit is designed for simple and rapid purification of total DNA from various sample sources including animal tissues and cells, blood, or gram negative bacteria. Purified DNA is free of contaminants and enzyme inhibitors and is highly suited for downstream applications such as PCR, blotting, and enzymatic reactions.

## Kit contents

Contents*			Storage conditions
Cat. No.	MG-P-005-50	MG-P-005-200	
No. of preparation	50	200	
Mini column and collecting tubes	50 ea	200 ea	Room Temperature
Buffer TL	15 mL	45 mL	
Buffer BL	15 mL	45 mL	
Buffer AW (conc.)**	18 mL	70 mL	
Buffer TW (conc.)**	8 mL	30 mL	
Buffer EA***	15 mL	50 mL	
Proteinase K solution (20 mg/mL)****	1.05 mL	4.2 mL	4°C

\* All components of this kit except Proteinase K solution should be stored at room temperature (15–25°C). Long exposure to heat source can deteriorate the performance of kit significantly.

\*\* Wash buffers are provided as concentrates. Ethanol must be added before first use as the indication on the bottle labels.

\*\*\* 10 mM Tris-HCl, pH 9.0, 0.5 mM EDTA

\*\*\*\* Proteinase K should be stored under 4°C on arrival for conservation of activity. Delivery time would not lead to noticeable loss of the enzyme's activity. Proteinase K solution can be stored at 4°C for a year without significant decrease of performance.

## Product use limitations

Axen™ Total DNA mini Kit is intended for research uses only. This kit is not intended for diagnosis or treatment for human. All due care and attention should be exercised in the handling of the products.

## Safety information

Axen™ Total DNA mini Kit contains irritants which are harmful when in contact with skin or eyes, or inhaled or swallowed. Care should be taken when handling this product. Always wear gloves and eye protection, and follow standard safety precautions.

Buffer BL and buffer AW contain chaotropes which can form highly reactive compounds when combined with bleach. Do NOT add bleach or acidic solutions directly to the sample-preparation waste.

## Preparation of ethanol-added buffer

Buffer AW and TW are provided as concentrate. Absolute ethanol (ACS grade or better) should be added as below before first use.

Buffer name	Included in	Volume of contents	Ethanol to be added	The final volume
AW	MG-P-005-50	18 mL	18 mL	36 mL
	MG-P-005-200	70 mL	70 mL	140 mL
TW	MG-P-005-50	8 mL	32 mL	40 mL
	MG-P-005-200	30 mL	120 mL	150 mL

## Protocol for preparation of total DNA

### *To check before start*

- Prepare water bath or heat block at 56°C
- Required materials; Absolute ethanol  
1.5 mL or 2 mL centrifuge tube,  
50 mL conical tube (Saliva and mouthwash)  
1x PBS (Swab, Saliva and mouthwash)
- Required optional materials; 100 mg/mL RNase A solution

**!!! Before first use, add the indicated volume of absolute ethanol to buffer AW and TW. (Refer to the *previous page*)**

### ● *Sample Preparation Procedure*

#### A. Cultured Cells

1. Harvest up to  $5 \times 10^6$  cultured cells in a 1.5 mL tube by centrifugation.
2. Discard the supernatant as much as possible.
3. Apply 200 uL of buffer TL and resuspend the cell pellet completely.
  - *For efficient resuspending in buffer TL, it is helpful to loosen and disrupt the pellets by flickering or vortexing before the addition of buffer TL.*
4. Pipet 20 uL of Proteinase K solution (20 mg/mL) into the tube and mix completely by vortexing or pipetting.
5. Go to 'DNA Purification procedure'

## **B. Animal Tissues**

- 1. Place up to 20 mg of a ground or minced tissue into a 1.5 mL tube.**
  - *Grind the tissue to a fine powder with liquid nitrogen in a pre-chilled mortar and pestle, or the tissue sample can be applied as minced with sterile sharp blade or scalpel as small as possible.*
  - *Alternatively, a soft tissue can be homogenized in buffer TL using a micro-homogenizer.*
  - *Do not exceed 10 mg for spleen or kidney.*
- 2. Apply 200  $\mu$ L of buffer TL and 20  $\mu$ L of Proteinase K solution (20 mg/mL) into the tube.**
- 3. Mix completely by vortexing or pipetting.**
  - *It is essential for good result to mix the components thoroughly to make a homogenate.*
- 4. Incubate the tube at 56°C until the tissue is completely lysed.**
  - *Complete lysis time varies from 10 mins to 3 hrs depending on the type of tissue. The lysate usually becomes translucent without any particulates after complete lysis. Overnight lysis does not influence the preparation.*
  - *For efficient lysis, do vortex occasionally during incubation. Use of a specialized instrument such as thermo-mixer will accelerate the lysis.*
  - *Make sure the tissue is completely digested.*
- 5. Go to 'DNA Purification procedure'**

### **C. Whole Blood (Mammalian)**

- 1. Place up to 200  $\mu$ L of whole blood sample into a 1.5 mL tube.**
  - *If a less sample than 200  $\mu$ L is used, adjust the sample volume to 200 $\mu$ L with buffer TL.*
- 2. Pipet 20  $\mu$ L of Proteinase K solution (20 mg/mL) into the tube and mix completely by vortexing or pipetting.**
- 3. Go to 'DNA Purification procedure'**

### **D. Nucleated Blood**

- 1. Place up to 10  $\mu$ L of nucleated blood sample with 190  $\mu$ L of buffer TL into a 1.5 mL tube.**
- 2. Pipet 20  $\mu$ L of Proteinase K solution (20 mg/mL) into the tube and mix completely by vortexing or pipetting.**
- 3. Go to 'DNA Purification procedure'**

### **E. Gram Negative Bacteria**

- 1. Harvest up to  $2 \times 10^9$  bacterial cells in a 1.5 mL tube by centrifugation.**
- 2. Discard the supernatant as much as possible.**
- 3. Apply 200  $\mu$ L of buffer TL and resuspend the pellet completely.**
  - *For efficient resuspending in buffer TL, it is helpful to loosen and disrupt the pellets by flickering or vortexing before the addition of buffer TL.*
- 4. Pipet 20  $\mu$ L of Proteinase K solution (20 mg/mL) into the tube and mix completely by vortexing or pipetting.**
- 5. Go to 'DNA Purification procedure'**

## **F. Buccal Swab**

- ☞ Buffer BL for the shortfalls due to the need of additional buffer BL for this procedure can be purchased separately.

### **1. Scrape the swab formly more than 5~6 times against the inside of cheek.**

- *To avoid contamination from other materials, ensure that the person who provides the sample has not taken any food or drink in 30 mins prior to sample collecting.*

### **2. Place the swab in 2.0 mL microcentrifuge tube. Clip off handle of brush with sterile sharp blade or wire cutter. Add 400 uL of 1X PBS to the tube.**

### **3. (*Optional*) If RNA-free DNA is required, add 4 uL of RNase A solution (100 mg/mL, not provided), vortex to mix, and incubate for 2 mins at room temperature.**

### **4. Apply 20 uL of Proteinase K solution (20 mg/mL, provided) and 400 uL of buffer BL to the sample. Vortex vigorously to mix completely.**

### **5. Incubate at 56°C for 10 mins.**

- *Briefly centrifuge to remove any drops from inside the lid.*

### **6. Add 400 uL of absolute ethanol (not provided) to the lysate and mix well by vortexing.**

- *Briefly centrifuge to remove any drops from inside the lid.*

### **7. Go to step 5 of 'DNA Purification procedure'**

## **G. Saliva and Mouthwash**

- 1. Collect 10 mL of mouthwash in a 50 mL conical tube or collect 1 mL of saliva by spitting in a 50 mL conical tube. If saliva is used, add 5 mL of 1x PBS to the sample and vortex to mix.**
  - *To avoid contamination from other materials, ensure that the person who provides the sample has not taken any food or drink in 30 mins prior to sample collecting.*
- 2. Centrifuge at 2,000 x g (3,000 rpm) for 5 mins to pellet cells. Immediately and carefully decant the supernatant to prevent loose cell pellets. Resuspend completely the pellets in 200 uL of buffer TL.**
  - *If the pellets are loose, repeat centrifugation.*
- 3. (*Optional*) If RNA should be removed, add 4 uL of RNase A solution (100 mg/mL, not provided) to the sample tube, vortex to mix and incubate for 2 mins at room temperature.**
- 4. Apply 20 uL of Proteinase K solution (20 mg/mL, provided) into the tube.**
- 5. Go to step 2 of 'DNA Purification procedure'**

## H. Hair

### 1. Prepare 'Buffer A' as follows:

#### **Buffer A;**

10 mM Tris-HCl, pH 8.0  
10 mM EDTA  
100 mM NaCl  
2% SDS  
40 mM DTT

- *Add DTT immediately before use, because it oxidizes quickly in a aqueous solutions.*

### 2. Collect hair sample in a 1.5 mL microcentrifuge tube.

- *The amount of starting material should not exceed 25 mg. It is recommended to use 0.5 ~ 1 cm from the root ends of plucked hair samples.*

### 3. Add 180 uL of prepared buffer A and 20 uL of Proteinase K solution (20 mg/mL, provided) to the tube, and vortex to mix thoroughly.

### 4. Incubate at 56°C for at least 1 hour until the sample is dissolved. Spin down briefly to remove any drops from inside of the lid.

- *Invert the tube occasionally to disperse the sample, or place on a rocking platform.*
- *Hair follicles should be completely dissolved, however hair shaft may be not dissolved completely and this residual solid materials will not affect DNA recovery.*

### 5. Go to 'DNA Purification procedure'



## I. Sperm

### 1. Prepare 'Buffer A2x' as follows:

#### **Buffer A2x;**

20 mM Tris-HCl, pH 8.0

20 mM EDTA

200 mM NaCl

4% SDS

80 mM DTT

- *Add DTT immediately before use, because it oxidizes quickly in a aqueous solutions.*
- 2. **Place 100 uL of sperm in a 1.5 mL microcentrifuge tube. Add 100 uL of buffer A2x and 20 uL of Proteinase K solution (20 mg/mL, provided) to the tube. Mix thoroughly by vortexing.**
- 3. **Incubate at 56°C until the sample is dissolved completely. Spin down briefly to remove any drops from inside of the lid.**
- 4. **Go to 'DNA Purification procedure'**

## ● **DNA Purification Procedure**

1. (**Optional**) If RNA should be removed, add 4 uL of RNase A solution (100 mg/mL, not provided) to the sample tube, vortex to mix and incubate for 2 mins at room temperature.
2. Add 200 uL of buffer BL and vortex to mix thoroughly.
  - *It is essential for good result to mix the components thoroughly to make a homogenate.*
3. Incubate the tube at 56°C for 10 mins.
4. Apply 200 uL of absolute ethanol and briefly vortex to mix.
5. Transfer up to 700 uL of the mixture into a spin column, centrifuge at 10,000 xg for 1 min, empty the collecting tube and re-insert the spin column back to the tube.
6. Repeat step 5 until all the remaining mixture has been applied.
7. Apply 600 uL of buffer AW to the spin column, centrifuge at 10,000 xg for 30 secs, empty the collecting tube and re-insert the spin column back to the tube.
8. Apply 600 uL of buffer TW to the spin column, centrifuge at 10,000 xg for 30 secs, empty the collecting tube and re-insert the spin column back to the tube.
9. Centrifuge the tube for 1 min at full speed for drying the membrane and transfer the spin column into a new 1.5 mL tube.
  - *Residual ethanol in eluate may interfere with the subsequent reactions. If carryover of wash buffer occurs, centrifuge again for 1 min before proceeding to next step.*
10. Apply 50 ~ 200 uL of buffer EA to the center of spin column membrane and incubate for 1 min at room temperature.
  - *Make sure that the buffer EA is dispensed directly onto the center of column membrane.*
  - *Elution volume can be decreased to 50 uL for higher concentration of the eluted DNA, but this will slightly decrease in overall DNA yield. If maximum recovery of DNA is preferred or the starting materials contain large amount of DNA, the elution volume can be increased to 200 uL.*
11. Centrifuge at full speed (>13,000 xg) for 1 min to elute the DNA.

## Trouble shooting guide

Facts	Possible causes	Recommendations
Low or no recovery	<i>Low cells in the starting sample</i>	Some samples may have low population of cells, and some whole blood may contain low numbers of white blood cells. Increase the sample volumes and load the binding mixture into the spin column several times. Reduce the elution volume to minimum. When the cell mass is very low, it is recommended to use a carrier RNA, which can be used by dissolving at 5~10 ug/mL in buffer BL.
	<i>Too many cells in the starting sample</i>	The sample amount over the maximum capacity will lead to poor lysis of cells, resulting in significantly low recovery. Reduce the amount of starting sample or increase the volume of buffers proportionally.
	<i>Inefficient lysis</i>	Some hard tissues such as heart and muscle take a long time to be completely lysed, and well shredded tissue will shorten this time. At start, pulverize the tissue to get homogenate. Using the thermo-mixer can greatly accelerate the lysis. Completely lysed mixture will not have any particulate on it. Degenerated Proteinase K may also cause inefficient lysis. Proteinase K should be stored under 4°C for maintenance of proper activity.
	<i>Improper eluent</i>	An eluent other than buffer EA can be employed as user's needs. The optimal eluting condition should be under low salt concentration with alkaline pH (7.0<pH<9.0). When water or other buffer was used as eluent, ensure that condition.
	<i>Ethanol was not added to the buffer</i>	Buffers AW and TW are provided as concentrate. Ethanol must be added to these buffers before first use. Otherwise, the result will be significantly poor.

Low purity or the colored residue on spin column membrane	<i>Inefficient lysis</i>	Inefficient lysis of sample may result in low DNA purity, and it may be due to the reasons described upper section, ' <i>Inefficient lysis</i> '.
	<i>Incomplete removal of hemoglobin</i>	In case of DNA preparation from a certain animal blood or from an improperly stored blood, it is hard to remove hemoglobin from their blood. Carry out additional wash step with buffer AW before washing with buffer TW.
High $A_{260}/A_{280}$ ratio	<i>RNA contamination</i>	RNA can inhibit some downstream enzymatic reactions, but not PCR itself. If RNA-free DNA is required, carry out RNase A treatment on step 1 of 'DNA Purification Procedure'.
Clogging of spin column	<i>Inefficient lysis</i>	Small particulates in incomplete lysate can clog the pore of the filter membrane. Please refer to the recommendation for inefficient lysis above.
Degraded DNA	<i>Starting sample is too old or improperly stored</i>	Too old or improperly stored samples have degraded DNA. Use fresh sample.
Enzymatic reaction is not performed well with purified DNA	<i>High salt concentration in eluate</i>	Ensure that wash step is carried out just in accordance with the protocols. Additional wash step may help to remove high salt in eluate.
	<i>Low purity of DNA eluate</i>	Please refer to the recommendation for 'low purity' above.
	<i>Residual ethanol in eluate</i>	Ensure that the wash step on the protocol is performed properly. The spin column membrane should be dried completely before eluting DNA. Perform an additional centrifugation to dry the membrane, if needed.