


Axen™ Taq DNA Polymerase



Axen™ Taq DNA Polymerase

Storage 1year at -20℃ | For research use only

1. Features

- 5'→3' exonuclease activity : Yes
- 3'→5' exonuclease activity : No
- Amplification size : < 3 kb
- PCR products : 3'-dA overhang
- Works with GC rich targets.

2. Kit components

	MG-E-001-500	MG-1-001-1000
Taq DNA Polymerase (5 Unit/ μ L)	500U	1,000U
10X Taq Buffer ^[1]	1.2mL	2 x 1.2mL
5X Q GC Enhancer ^[2]	2mL	2 x 2mL
dNTP Mix, 10mM each	0.3mL	0.5mL
Water, nuclease free	2mL	2mL

^[1] Includes 20mM MgCl₂


^[2] For amplicons with > 65% GC content

3. Protocol

1) Reaction Mixture (For 20 μ L or 50 μ L reaction)

	20 μ L rxn	50 μ L rxn
10X Taq Buffer	2 μ L	5 μ L
5X Q GC Enhancer	4 μ L	10 μ L
dNTP Mix, 10mM each	0.4 μ L	1 μ L
Forward primer	0.1–1.0 μ M	
Reverse primer	0.1–1.0 μ M	
Template DNA	5–100ng genomic DNA 0.01–10ng plasmid	
Taq DNA Polymerase	0.2 μ L (1U)	0.5 μ L (2.5U)
Water, nuclease free	add to 20 μ L	add to 50 μ L

2) Thermal cycling conditions



^[3] Optimal annealing temperature depends on the melting temperature of the primers and on the system used.

^[4] The recommended extension step is 1 minute for PCR products up to 2 kb. For longer products, the extension time should be prolonged by 1 minute/kb.